

## Instructions for use

# DHT ELISA

**REF****AA E-1700R****RUO**

For research  
use only –  
Not for use  
in diagnostic  
procedures

## **Dihydrotestosterone (DHT) ELISA**

### ***INTENDED USE***

For the direct quantitative determination of Dihydrotestosterone by an enzyme immunoassay in human serum.

### ***PRINCIPLE OF THE TEST***

The principle of the following enzyme immunoassay test follows the typical competitive binding scenario. Competition occurs between an unlabelled antigen (present in standards, controls and samples) and an enzyme-labelled antigen (conjugate) for a limited number of antibody binding sites on the microplate. The washing and decanting procedures remove unbound materials. After the washing step, the enzyme substrate is added. The enzymatic reaction is terminated by addition of the stopping solution. The absorbance is measured on a microtiter plate reader. The intensity of the colour formed is inversely proportional to the concentration of DHT in the sample. A set of standards is used to plot a standard curve from which the amount of DHT in samples and controls can be directly read.

### ***PROCEDURAL CAUTIONS AND WARNINGS***

1. Users should have a thorough understanding of this protocol for the successful use of this kit. Reliable performance will only be attained by strict and careful adherence to the instructions provided.
2. Control materials or serum pools should be included in every run at a high and low level for assessing the reliability of results.
3. When the use of water is specified for dilution or reconstitution, use deionized or distilled water.
4. In order to reduce exposure to potentially harmful substances, gloves should be worn when handling kit reagents and human specimens.
5. All kit reagents and specimens should be brought to room temperature and mixed gently but thoroughly before use. Avoid repeated freezing and thawing of reagents and specimens.
6. A standard curve must be established for every run.
7. The controls should be included in every run and fall within established confidence limits.
8. Improper procedural techniques, imprecise pipetting, incomplete washing as well as improper reagent storage may be indicated when assay values for the controls do not reflect established ranges.
9. When reading the microplate, the presence of bubbles in the wells will affect the optical densities (ODs). Carefully remove any bubbles before performing the reading step.
10. The substrate solution (TMB) is sensitive to light and should remain colourless if properly stored. Instability or contamination may be indicated by the development of a blue colour, in which case it should not be used.
11. When dispensing the substrate and stopping solution, do not use pipettes in which these liquids will come into contact with any metal parts.
12. To prevent contamination of reagents, use a new disposable pipette tip for dispensing each reagent, sample, standard and control.
13. Do not mix various lot numbers of kit components within a test and do not use any component beyond the expiration date printed on the label.
14. Kit reagents must be regarded as hazardous waste and disposed of according to national regulations.

### ***LIMITATIONS***

1. All the reagents within the kit are calibrated for the direct determination of DHT in human serum. The kit is not calibrated for the determination of DHT in saliva, plasma or other specimens of human or animal origin.
2. Do not use grossly hemolyzed, grossly lipemic, icteric or improperly stored serum.
3. Any samples or control sera containing azide or thimerosal are not compatible with this kit, as they may lead to false results.
4. Only standard A may be used to dilute any high serum samples. The use of any other reagent may lead to false results.

### ***SAFETY CAUTIONS AND WARNINGS***

#### ***POTENTIAL BIOHAZARDOUS MATERIAL***

Human serum that may be used in the preparation of the standards and controls has been tested and found to be nonreactive for Hepatitis B surface antigen and has also been tested for the presence of antibodies to HCV and Human Immunodeficiency Virus (HIV) and found to be negative. No test method however, can offer complete assurance that HIV, HCV and Hepatitis B virus or any infectious agents are absent.

The reagents should be considered a potential biohazard and handled with the same precautions as applied to any blood specimen.

#### ***CHEMICAL HAZARDS***

Avoid contact with reagents containing TMB, hydrogen peroxide and sulfuric acid. If contacted with any of these reagents, wash with plenty of water. TMB is a suspected carcinogen.

#### ***SPECIMEN COLLECTION AND STORAGE***

Approximately 0.2 ml of serum is required per duplicate determination. Collect 4–5 ml of blood into an appropriately labelled tube and allow it to clot. Centrifuge and carefully remove the serum layer. Store at 4°C

for up to 24 hours or at -10°C or lower if the analyses are to be done at a later date. Consider all human specimens as possible biohazardous materials and take appropriate precautions when handling.


### ***SPECIMEN PRETREATMENT***

This assay is a direct system; no specimen pretreatment is necessary.

### ***REAGENTS AND EQUIPMENT NEEDED BUT NOT PROVIDED***

1. Precision pipettes to dispense 50, 100, 150 and 300 µl
2. Disposable pipette tips
3. Distilled or deionized water
4. Microplate reader with a filter set at 450 nm and an upper OD limit of 3.0 or greater\*  
(see assay procedure step 10)

### ***REAGENTS PROVIDED***

- 1. AA E-0030** WASH-CONC 10x **Wash Buffer Concentrate** – Requires Preparation X10  
Contents: One bottle containing buffer with a non-ionic detergent and a non-mercury preservative.  
Volume: 50 ml/bottle  
Storage: Refrigerate at 2 - 8 °C  
Stability: 12 months or as indicated on label.  
Preparation: Dilute 1:10 in distilled or deionized water before use. If the whole plate is to be used dilute 50 ml of the wash buffer concentrate in 450 ml of water.
- 2. AA E-0055** SUBSTRATE **TMB Substrate** - Ready To Use.  
Contents: One bottle containing tetramethylbenzidine and hydrogen peroxide in a non-DMF or DMSO containing buffer.  
Volume: 16 ml/bottle  
Storage: Refrigerate at 2 - 8 °C  
Stability: 12 months or as indicated on label.
- 3. AA E-0080** STOP-SOLN **Stopping Solution** - Ready To Use.  
Contents: One bottle containing 1M sulfuric acid.  
Volume: 6 ml/bottle  
Storage: Refrigerate at 2 - 8 °C  
Stability: 12 months or as indicated on label.  
Hazards identification:   
H290 May be corrosive to metals.  
H314 Causes severe skin burns and eye damage.

#### 4. Standards and Controls- Ready To Use.

Listed below are approximate concentrations, please refer to vial labels for exact concentrations:

Cat. no.	Symbol	Standard	Concentration	Volume/Vial
AA E-1701	STANDARD A	Standard A	0 pg/ml	2.0 ml
AA E-1702	STANDARD B	Standard B	25 pg/ml	0.6 ml
AA E-1703	STANDARD C	Standard C	100 pg/ml	0.6 ml
AA E-1704	STANDARD D	Standard D	500 pg/ml	0.6 ml
AA E-1705	STANDARD E	Standard E	1000 pg/ml	0.6 ml
AA E-1706	STANDARD F	Standard F	2500 pg/ml	0.6 ml
AA E-1751	CONTROL 1	Control 1	Refer to vial labels for expected value and acceptable range!	0.6 ml
AA E-1752	CONTROL 2	Control 2		0.6 ml

Contents: DHT in a protein-based buffer with a non-mercury preservative. Prepared by spiking buffer with a defined quantity of DHT.

Storage: Refrigerate at 2 - 8 °C

Stability: 12 months in unopened vial or as indicated on label. Once opened, the standards and controls should be used within 14 days or aliquoted and stored frozen.

Avoid multiple freezing and thawing cycles.

#### 5. AA E-1713 ASSAY-BUFF Assay Buffer - Ready To Use.

Contents: One bottle containing a protein-based buffer with a non-mercury preservative.

Volume: 15 ml/vial

Storage: Refrigerate at 2 - 8 °C

Stability: 12 months or as indicated on label.

#### 6. AA E-1731 96 Rabbit Anti-DHT Antibody-Coated Break-Apart Well Microplate - Ready To Use.

Contents: One 96 well (12x8) polyclonal antibody-coated microplate in a resealable pouch with desiccant.

Storage: Refrigerate at 2 - 8 °C

Stability: 12 months or as indicated on label.

#### 7. AA E-1740 CONJUGATE-CONC 100x Dihydrotestosterone-Horseradish Peroxidase (HRP) Conjugate Concentrate – Requires Preparation X100

Contents: DHT-HRP conjugate in a protein-based buffer with a non-mercury preservative.

Volume: 200 µl/vial

Storage: Refrigerate at 2 - 8 °C

Stability: 12 months or as indicated on label.


Preparation: Dilute 1:100 in assay buffer before use (eg. 20 µl of conjugate concentrate in 2 ml of assay buffer). If the whole plate is to be used dilute 120 µl of conjugate concentrate in 12 ml of assay buffer.

Discard any that is left over.

## ASSAY PROCEDURE

Specimen Pretreatment: **None**.

All reagents must reach room temperature before use. Standards, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.

<b>1.</b>	Prepare <b>working solutions</b> of the <b>DHT-HRP conjugate</b> and <b>wash buffer</b> .
<b>2.</b>	Remove the required number of well strips. Reseal the bag and return any unused strips to the refrigerator.
<b>3.</b>	<b>Pipette 50 µl</b> of each <b>standard, control and specimen sample</b> into correspondingly labelled wells in duplicate.
<b>4.</b>	Pipette <b>100 µl</b> of the <b>conjugate working solution</b> into each well. <i>(We recommend using a multichannel pipette).</i>
<b>5.</b>	Gently shake the plate for <b>10 seconds</b> and incubate for <b>1 hour</b> at <b>room temperature</b> <i>(no shaking)</i> .
<b>6.</b>	Wash the wells <b>3 times</b> with <b>300 µl</b> of diluted wash buffer per well and tap the plate firmly against absorbent paper to ensure that it is dry <i>(The use of a washer is recommended)</i> .
<b>7.</b>	Pipette <b>150 µl</b> of <b>TMB substrate</b> into each well at timed intervals.
<b>8.</b>	Gently shake the plate for <b>10 seconds</b> and incubate for <b>10 – 15 minutes</b> at <b>room temperature</b> <i>(no shaking)</i> or until Standard A attains dark blue colour for desired OD.
<b>9.</b>	Pipette <b>50 µl</b> of <b>stopping solution</b> into each well at the same timed intervals as in step 7.
<b>10.</b>	Read the plate on a microwell plate reader at <b>450 nm</b> within 20 minutes after addition of the stopping solution.
	<i>If the OD exceeds the upper limit of detection or if a 450 nm filter is unavailable, a 405 or 415nm filter may be substituted. The optical densities will be lower, however, this will not affect the results of control samples.</i>

## CALCULATIONS


1. Calculate the mean optical density of each standard duplicate.
2. Draw a standard curve on semi-log paper with the mean optical densities on the Y-axis and the standard concentrations on the X-axis. If immunoassay software is being used, a 4-parameter or 5-parameter curve is recommended.
3. Calculate the mean optical density of each unknown duplicate.
4. Read the values of the unknowns directly off the standard curve.
5. If a sample reads more than 2500 pg/ml then dilute it with standard A at a dilution of no more than 1:8. The result obtained should be multiplied by the dilution factor.

## TYPICAL TABULATED DATA:

Standard	OD 1	OD 2	Mean OD	Value (pg/ml)
A	2.320	2.279	2.300	0
B	1.976	1.928	1.952	25
C	1.058	1.077	1.068	100
D	0.359	0.354	0.357	500
E	0.222	0.205	0.214	1000
F	0.131	0.128	0.130	2500
Unknown	0.515	0.507	0.511	300

**TYPICAL STANDARD CURVE**

Sample curve only. **Do not** use to calculate results.



**PERFORMANCE CHARACTERISTICS**

**SENSITIVITY**

The lower detection limit is calculated from the standard curve by determining the resulting concentration of the mean OD of standard A (based on 10 replicate analyses) minus 2 SD. Therefore, the sensitivity of the Direct Dihydrotestosterone ELISA kit is **6.0 pg/ml**.

**SPECIFICITY (CROSS REACTIVITY)**

The following compounds were tested for cross-reactivity with the Direct Dihydrotestosterone ELISA kit with dihydrotestosterone cross-reacting at 100%.

Steroid	% Cross Reactivity
Dihydrotestosterone	100
Testosterone	8.7
5 $\beta$ Dihydrotestosterone	2.0
Androstenedione	0.2

The following steroids were tested but cross-reacted at less than 0.01%: Dehydroepiandrosterone Sulfate, 17 $\beta$ -Estradiol, Estriol, Estrone, Progesterone, 17-OH Progesterone, Cortisol, and Pregnenolone.

**INTRA-ASSAY PRECISION**

Three samples were assayed ten times each on the same standard curve. The results (in pg/ml) are tabulated below:

Sample	Mean	SD	CV%
1	236.74	26.89	11.4
2	869.03	47.41	5.46
3	1008.14	39.36	3.90

**INTER-ASSAY PRECISION**

Three samples were assayed ten times over a period of four weeks. The results (in pg/ml) are tabulated below:

Sample	Mean	SD	CV%
1	280.88	34.07	12.1
2	721.40	54.20	7.5
3	1025.41	60.45	5.9

**RECOVERY**

Spiked samples were prepared by adding defined amounts of DHT to three serum samples. The results (in pg/ml) are tabulated below:

Sample	Obs.Result	Exp.Result	Recovery%
1 Unspiked	290.54	-	-
+117.53	361.51	408.07	88.6
+235.06	501.66	525.60	95.4
+470.13	744.81	760.67	97.9
2 Unspiked	324.75	-	-
+117.53	389.43	442.29	88.0
+235.06	505.23	559.81	90.3
+470.13	712.44	794.88	89.6
3 Unspiked	720.11	-	-
+117.53	758.13	837.64	90.5
+235.06	856.46	955.17	89.7
+470.13	1013.61	1190.24	85.1

**LINEARITY**

Three serum samples were diluted with Standard A. The results (in pg/ml) are tabulated below:

Sample	Obs.Result	Exp.Result	Recovery%
1	340.67	-	-
1:2	165.35	170.34	97.1
1:4	95.39	85.17	112.0
1:8	48.47	42.58	113.8
2	1086.01	-	-
1:2	508.58	543.00	93.7
1:4	232.11	271.50	85.5
1:8	114.95	135.75	84.7
3	1313.21	-	-
1:2	612.98	656.61	93.4
1:4	318.63	328.30	97.1
1:8	134.98	164.15	82.2

**COMPARATIVE STUDIES**

The Direct Dihydrotestosterone ELISA kit (Kit A) was compared with a competitors coated tube RIA kit (Kit B). The results (in pg/ml) are tabulated below:

Group	N	Kit A Mean	Kit B Mean
Females	10	95.5	99.0
Males	10	280.0	252.0

**EXPECTED NORMAL VALUES**







As for all clinical assays each laboratory should collect data and establish their own range of expected normal values.

Group	Range (pg/ml)
Females:	
Pre-menopausal	24-368
Postmenopausal	10-181
Males	250-990

**REFERENCES**

1. Bassett RM. A Simple Chromatographic Method for the Radioimmunoassay of Four Androgenic Steroids. *Med Lab Sci.* 1980; 37(1):31-8.
2. Baxendale PM, et al. Plasma and Salivary Androstenedione and Dihydrotestosterone in Women with Hyperandrogenism. *Clin Endocrinol (Oxf).* 1983; 18(5):447-57.
3. Brooks RV, Androgens. Physiology and Pathology. In: Makin HL J, ed., *Biochemistry of Steroid Hormones*, 2nd ed., Oxford Blackwell Scientific Publications, 1984; 565.
4. Cameron EHD, et al., eds. *Steroid Immunoassay: Proceedings of the Fifth Tenovus Workshop*, Cardiff, April 1974. Cardiff: Alpha Omega Publishing; 1975.
5. Dunn JF, et al. Transport of Steroid Hormones: Binding of 21 Endogenous Steroids to Testosterone-Binding Globulin and Corticosteroid-Binding Globulin in Human Plasma. *J Clin Endocrinol Metab.* 1981; 53(1):58-68.
6. Hammond GL, et al. The Simultaneous Radioimmunoassay of Seven Steroids in Human Spermatic and Peripheral Venous Blood. *J Clin Endocrinol Metab.* 1977; 45(1):16-24.
7. Ito T, Horton R. Dihydrotestosterone in Human Peripheral Plasma. *J Clin Endocrinol Metab.* 1970; 31(4):362-8.
8. Mean F, et al. Study of the Binding of Dihydrotestosterone, Testosterone and Oestradiol with Sex Hormone Binding Globulin. *Clin Chim Acta.* 1977; 80(1):171-80.
9. Mooradian AD, et al. The Biological Actions of Androgens. *Endocr Rev.* 1987; 8(1):1-28.
10. Pazzagli M, et al. Radioimmunoassay of Plasma Dihydrotestosterone in Normal and Hypogonadal Men. *Clin Endocrinol (Oxf).* 1975; 4(5):513-20.
11. Wakelin K, et al. Relationship of 5 $\beta$ dihydrotestosterone and 5 $\alpha$ dihydrotestosterone to testosterone in health and disease. *J Endocrinol.* 1980; 87:450.
12. Wang C, et al. Salivary Androgens in Hirsutism: Are They of Use in Routine Evaluation? *Ann Clin Biochem.* 1986; 23(Pt 5):590-5.
13. Kricka LJ. Human Anti-Animal Antibody Interferences in Immunological Assays. *Clin Chem.* 1999; 45(7):942-56.
14. Check JH, et al. Falsely Elevated Steroidal Assay Levels Related to Heterophile Antibodies Against Various Animal Species. *Gynecol Obstet Invest.* 1995; 40(2):139-40.

**Symbols:**

	Storage temperature		Manufacturer		Contains sufficient for <n> tests
	Expiry date	<b>LOT</b>	Batch code		
	Consult instructions for use	<b>CONT</b>	Content		
	Caution	<b>REF</b>	Catalogue number	<b>RUO</b>	For research use only!