**IMMUNOASSAYS AND SERVICES**

**BIOGENIC AMINES & NEUROSCIENCE | ENDOCRINOLOGY | FOOD SAFETY**

LABOR DIAGNOSTIKA NORD GmbH & Co.KG | Am Eichenhain 1 | 48531 Nordhorn | Germany | Tel. +49 5921 8197-0 | Fax +49 5921 8197-222 | info@ldn.de | www.ldn.de

---

**Instructions for use**  
**Histamine ELISA** **Fast Track**

**REF**

**BA E-1700**

**IVD**

## **1. Introduction**

### **1.1 Intended use and principle of the test**

Enzyme Immunoassay for the quantitative determination of Histamine in urine and whole blood.

For research use only (RUO) this kit can also be used for the quantitative determination of Histamine cell culture samples. Please refer to appendix 1.

First, Histamine is quantitatively acylated. The subsequent competitive ELISA kit uses the microtiter plate format. The antigen is bound to the solid phase of the microtiter plate. The acylated analyte concentrations in the standards, controls and samples and the solid phase bound analyte compete for a fixed number of antibody binding sites. After the system is in equilibrium, free antigen and free antigen-antibody complexes are removed by washing. The antibody bound to the solid phase is detected by an anti-goat IgG-peroxidase conjugate using TMB as a substrate. The reaction is monitored at 450 nm.

Quantification of unknown samples is achieved by comparing their absorbance with a standard curve prepared with known standard concentrations.

### **1.2 Clinical application**

Histamine belongs to the biogenic amines and is synthesized by decarboxylation from the amino acid histidine. It is synthesized by mast cells, basophils, platelets, histaminergic neurons, and enterochromaffine cells, where it is stored intracellularly in vesicles and released on stimulation.

Histamine acts by binding to its 4 receptors (H1R, H2R, H3R and H4R) on target cells in various tissues. It causes smooth muscle cell contraction, vasodilatation, increased vascular permeability and mucus secretion, tachycardia, alterations of blood pressure, and arrhythmias.

In humans, histamine is one of the most important mediators and takes part in the initial phase of an anaphylactic reaction ("immediate type" allergy).

Of clinical interest is also the quantification of the histamine *release* from basophilic leucocytes in allergies.

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as under point "Procedural cautions, guidelines and warnings". Any laboratory result is only a part of the total clinical picture of the patient.

Only in cases where the laboratory results are in an acceptable agreement with the overall clinical picture of the patient it can be used for therapeutic consequences.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

## **2. Procedural cautions, guidelines, warnings and limitations**

### **2.1 Procedural cautions, guidelines and warnings**

- (1) This kit is intended for professional use only. Users should have a thorough understanding of this protocol for the successful use of this kit. Only the test instruction provided with the kit is valid and has to be used to run the assay. Reliable performance will only be attained by strict and careful adherence to the instructions provided.
- (2) This assay was validated for certain types of sample as indicated in *Intended Use* (please refer to Chapter 1). Any off-label use of this kit is in the responsibility of the user and the manufacturer cannot be held liable.
- (3) The principles of Good Laboratory Practice (GLP) have to be followed.
- (4) In order to reduce exposure to potentially harmful substances, wear lab coats, disposable protective gloves and protective glasses where necessary.
- (5) All kit reagents and specimens should be brought to room temperature and mixed gently but thoroughly before use. Avoid repeated freezing and thawing of reagents and specimens.
- (6) For dilution or reconstitution purposes, use deionized, distilled, or ultra-pure water.
- (7) The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch with desiccant and used in the frame provided.
- (8) Duplicate determination of sample is highly recommended to be able to identify potential pipetting errors.
- (9) Once the test has been started, all steps should be completed without interruption. Make sure that the required reagents, materials and devices are prepared ready at the appropriate time.
- (10) Incubation times do influence the results. All wells should be handled in the same order and time intervals.
- (11) To avoid cross-contamination of reagents, use new disposable pipette tips for dispensing each reagent, sample, standard and control.
- (12) A standard curve must be established for each run.
- (13) The controls should be included in each run and fall within established confidence limits. The confidence limits are listed in the QC-Report provided with the kit.
- (14) Do not mix kit components with different lot numbers within a test and do not use reagents beyond expiry date as shown on the kit labels.

- (15) Avoid contact with Stop Solution containing 0.25 M H<sub>2</sub>SO<sub>4</sub>. It may cause skin irritation and burns. In case of contact with eyes or skin, rinse off immediately with water.
- (16) TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them.
- (17) For information on hazardous substances included in the kit please refer to Material Safety Data Sheet (MSDS). The Material Safety Data Sheet for this product is made available directly on the website of the manufacturer or upon request.
- (18) The expected reference values reported in this test instruction are only indicative. It is recommended that each laboratory establishes its own reference intervals.
- (19) The results obtained with this test kit should not be taken as the sole reason for any therapeutic consequence but have to be correlated to other diagnostic tests and clinical observations.
- (20) Kit reagents must be regarded as hazardous waste and disposed of according to national regulations.

## 2.2 Limitations

Any inappropriate handling of samples or modification of this test might influence the results.

### 2.2.1 Interfering substances

#### 24-hour urine

Please note the sample preparation! If the percentage of the final concentration of acid is too high, the buffer capacity of the Acylation Buffer is insufficient. As a consequence histamine will not be acylated quantitatively.

### 2.2.2 Drug interferences

There are no known substances (drugs) which ingestion interferes with the measurement of histamine level in the sample.

### 2.2.3 High-Dose-Hook effect



No hook effect was observed in this test.

## 3. Storage and stability

Store the unopened reagents at 2 - 8 °C until expiration date. Do not use components beyond the expiry date indicated on the kit labels. Once opened the reagents are stable for 1 month when stored at 2 - 8 °C. Once the resealable pouch has been opened, care should be taken to close it tightly with desiccant again.

## 4. Materials

### 4.1 Contents of the kit

<b>BA D-0024</b>		<b>Reaction Plate</b> - Ready to use
Contents:	1 x 96 well plate, empty in a resealable pouch	
<b>BA E-0030</b>		<b>Wash Buffer Concentrate</b> - Concentrated 50x
Contents:	Buffer with a non-ionic detergent and physiological pH	
Volume:	1 x 20 ml/vial, light purple cap	
<b>BA E-1240</b>		<b>Enzyme Conjugate</b> - Ready to use
Contents:	Donkey anti-goat immunoglobulins conjugated with peroxidase	
Volume:	1 x 12 ml/vial, red cap	
<b>BA E-0055</b>		<b>Substrate</b> - Ready to use
Contents:	Chromogenic substrate containing tetramethylbenzidine, substrate buffer and hydrogen peroxide	
Volume:	1 x 12 ml/black vial, black cap	
<b>BA E-0080</b>		<b>Stop Solution</b> - Ready to use
Contents:	0.25 M sulfuric acid	
Volume:	1 x 12 ml/vial, light grey cap	
<b>BA E-1031</b>		<b>Histamine Microtiter Strips</b> - Ready to use
Contents:	1 x 96 well (12x8) antigen precoated microwell plate in a resealable foil pouch with desiccant	

**BA E-1210** **HIS-AS** **Histamine Antiserum** - Ready to use

Contents: Goat anti-histamine antibody, blue coloured

Volume: 1 x 12 ml/vial, blue cap

**Standards and Controls** - Ready to use

Cat. no.	Component	Colour/Cap	Concentration ng/ml	Concentration nmol/l	Volume/ Vial
<b>BA E-1001</b>	<b>STANDARD A</b>	white	0	0	4 ml
<b>BA E-1002</b>	<b>STANDARD B</b>	light yellow	0.5	4.5	4 ml
<b>BA E-1003</b>	<b>STANDARD C</b>	orange	1.5	13.5	4 ml
<b>BA E-1004</b>	<b>STANDARD D</b>	dark blue	5	45	4 ml
<b>BA E-1005</b>	<b>STANDARD E</b>	light grey	15	135	4 ml
<b>BA E-1006</b>	<b>STANDARD F</b>	black	50	450	4 ml
<b>BA E-1051</b>	<b>CONTROL 1</b>	light green	Refer to QC-Report for expected value and acceptable range!		4 ml
<b>BA E-1052</b>	<b>CONTROL 2</b>	dark red			4 ml

Conversion: Histamine (ng/ml) x 9 = Histamine (nmol/l)

Contents: Acidic buffer spiked with defined quantity of Histamine

**BA E-1211** **ACYL-BUFF** **Acylation Buffer** - Ready to use

Contents: TRIS-buffer

Volume: 2 x 12 ml/vial, brown cap

**BA E-1212** **ACYL-REAG** **Acylation Reagent** - Ready to use

Contents: Acylation reagent containing DMSO

Volume: 2 x 1.5 ml/vial

**BA E-0041** **DILUENT** **Diluent** - Ready to use

Contents: Acidic buffer with non-mercury preservatives

Volume: 1 x 22 ml/vial, white cap

**4.2 Additional materials and equipment required but not provided in the kit**

- Calibrated precision pipettes to dispense volumes between 20 - 200 µl
- Microtiter plate washing device (manual, semi-automated or automated)
- ELISA reader capable of reading absorbance at 450 nm and if possible 620 - 650 nm
- Absorbent material (paper towel)
- Water (deionized, distilled, or ultra-pure)
- Vortex mixer

⚠ The assay can be performed with or without the usage of a microtiter plate shaker. If a shaker is used, it should have the following characteristics: shaking amplitude 3 mm, capable of approx. 600 rpm.

⚠ For the determination of total histamine in whole blood, a **Release Buffer** is necessary!  
(Please ask your local supplier for ordering details)

**BA E-1726** **RELEASE-BUFF** **Release Buffer** - Ready to use

Contents: Buffer with physiological pH

Volume: 1 x 250 ml/vial, white cap

**5. Sample collection and storage****Urine**

Spontaneous urine or 24-hour urine, collected in a bottle containing 10 - 15 ml of 6 M HCl, may be used. *If 24-hour urine is being used please record the total volume of the collected urine.* If the percentage of the final concentration of acid is too high, the buffer capacity of the Diluent is insufficient. As a consequence interfering factors are not extracted quantitatively.

Storage: up to 6 hours at 2 - 8 °C; for longer periods (up to 6 months) at -20 °C.

Repeated freezing and thawing should be avoided.

Avoid exposure to direct sunlight.

### Whole Blood

Whole blood is collected into a tube (e.g. Monovette™ or Vacuette™) containing heparin as anti-coagulant. The samples can be stored for up to 24 hours at room temperature. Please do not keep the samples refrigerated, this will lead to clotting of the leucocytes. Avoid direct sunlight.

## 6. Test procedure

Allow all reagents and samples to reach room temperature and mix thoroughly by gentle inversion before use. Duplicate determinations are recommended.

The binding of the antisera and of the enzyme conjugate and the activity of the enzyme are temperature dependent, and the absorption values may vary if a thermostat is not used. The higher the temperature, the higher the absorption values will be. Varying incubation times will have similar influences on the absorbance. The optimal temperature during the Enzyme Immunoassay is between 20 – 25 °C.

### 6.1 Preparation of reagents

#### Wash Buffer

Dilute the 20 ml Wash Buffer Concentrate with water (deionized, distilled, or ultra-pure) to a final volume of 1000 ml.

Storage: 1 month at 2 – 8 °C

#### Acylation Reagent

The Acylation Reagent has a freezing point of 18.5 °C. To ensure that the Acylation Reagent is liquid when being used, it must be ensured that the Acylation Reagent has reached room temperature and forms a homogeneous, crystal-free solution.

#### Whole Blood

For the determination of total histamine in whole blood, dilute heparinized whole blood 1 + 20 with **Release Buffer** (BA E-1726) and incubate for **10 min** at **90 °C** (e.g. 50 µl whole blood plus 1 ml Release Buffer)

Cool down the samples for **10 min** at **2 – 8°C**

Centrifuge for 10 min at 700 x g (*the brake should be switched-off!*)

### 6.2 Sample preparation and acylation

1.	Pipette <b>100 µl</b> of <b>standards, controls</b> and <b>20 µl</b> of <b>urine samples</b> and <b>100 µl</b> of <b>perpetrated whole blood samples</b> ( <i>refer to 6.1.</i> ) into the respective wells of the <b>Reaction Plate</b> .
2.	Add <b>80 µl</b> of <b>Diluent</b> to the wells with <b>urine samples</b> .
3.	Add <b>25 µl</b> of <b>Acylation Reagent</b> ( <i>refer to 6.1</i> ) to all wells.
4.	Pipette <b>200 µl</b> of <b>Acylation Buffer</b> into all wells.
5.	Incubate <b>15 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm). <b>Alternatively without shaker: shake the Reaction Plate shortly by hand and incubate 15 min at RT (20 – 25 °C).</b>
	<b>Take 25 µl for the subsequent ELISA</b>

### 6.3 Histamine ELISA

1.	Pipette <b>25 µl</b> of the <b>acylated standards, controls and samples</b> into the wells of the <b>Histamine Microtiter Strips</b> .
2.	Pipette <b>100 µl</b> of the <b>Histamine Antiserum</b> into all wells.
3.	Incubate <b>30 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm). <i><b>Alternatively without shaker:</b> shake the <b>Histamine Microtiter Strips</b> shortly by hand and incubate for <b>40 min</b> at <b>RT</b> (20 – 25 °C).</i>
4.	Discard or aspirate the contents of the wells and wash the plate <b>3 x</b> by adding <b>300 µl</b> of <b>Wash Buffer</b> , <b>discarding</b> the content and <b>blotting dry each time</b> by tapping the inverted plate on absorbent material.
5.	Pipette <b>100 µl</b> of the <b>Enzyme Conjugate</b> into all wells.
6.	Incubate for <b>10 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm). <i><b>Alternatively without shaker:</b> incubate for <b>20 min</b> at <b>RT</b> (20 – 25 °C).</i>
7.	Discard or aspirate the contents of the wells and wash the plate <b>3 x</b> by adding <b>300 µl</b> of <b>Wash Buffer</b> , <b>discarding</b> the content and <b>blotting dry each time</b> by tapping the inverted plate on absorbent material.
8.	Pipette <b>100 µl</b> of the <b>Substrate</b> into all wells.
9.	Incubate for <b>15 ± 2 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm). <i><b>Alternatively without shaker:</b> incubate for <b>15 ± 2 min</b> at <b>RT</b> (20 – 25 °C).</i>
	<b>Avoid exposure to direct sunlight!</b>
10.	Add <b>100 µl</b> of the <b>Stop Solution</b> to each well and shake the microtiter plate shortly by hand to ensure a homogeneous distribution of the solution.
11.	<b>Read</b> the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to <b>450 nm</b> (if available a reference wavelength between 620 nm and 650 nm is recommended).

### 7. Calculation of results

Measuring range		<b>Histamine</b>
	Controls	0.12 - 50 ng/ml
	Urine	0.6 - 250 ng/ml
	Whole blood	2.52 - 1050 ng/ml

The standard curve is obtained by plotting the absorbance readings (calculate the mean absorbance) of the standards (linear, y-axis) against the corresponding standard concentrations (logarithmic, x-axis). Use a non-linear regression for curve fitting (e.g. spline, 4- parameter, akima).

 *This assay is a competitive assay. This means: the OD-values are decreasing with increasing concentrations of the analyte. OD-values found below the standard curve correspond to high concentrations of the analyte in the sample and have to be reported as being positive.*

#### Controls

The concentrations of the **controls** can be read directly from the standard curve.

#### Whole blood samples

The read concentrations of **histamine in whole blood** have to be **multiplied by 21**.

#### Urine samples

The read concentrations of **histamine in urine** have to be **multiplied by 5**.

The total amount of Histamine excreted in urine during 24 h is calculated as following:

$$\mu\text{g}/24\text{h} = \mu\text{g}/\text{l} \times \text{l}/24\text{h}$$

#### Conversion

$$\text{Histamine (ng/ml)} \times 9 = \text{Histamine (nmol/l)}$$

#### Expected reference values

It is strongly recommended that each laboratory should determine its own reference values.

<b>Whole blood (total Histamine)</b>	<b>24 hour-urine</b>	<b>Spontaneous urine</b>
20 - 200 ng/ml	< 45 µg/d	< 45 µg/g creatinine


## 7.1 Quality control

It is recommended to use control samples according to national regulations. Use controls at both normal and pathological levels. The kit or other commercially available controls should fall within established confidence limits. The confidence limits of the kit controls are listed in the QC-Report.

## 7.2 Typical standard curve

Example, do not use for calculation!

## 8. Assay characteristics

Analytical Specificity (Cross Reactivity)	Substance	Cross Reactivity (%)
		Histamine
	3-Methyl-Histamine	0.1
	Tyramine	0.01
	L-Phenylalanine	< 0.001
	L-Histidine	< 0.001
	L-Tyrosine	< 0.001
	Tryptamine	< 0.001
	5-Hydroxy-Indole-Acetic Acid	< 0.001
	Serotonine	< 0.001

<b>Analytical Sensitivity</b> (Limit of Detection)	Histamine	0.2 ng/ml	Mean signal (Zero-Standard) + 2SD
---	-----------	-----------	-----------------------------------

Precision					
Inter-Assay Variation, n = 13			Intra-Assay Variation, n = 39		
Sample	Mean ± SD (ng/ml)	CV (%)	Sample	Mean ± SD (ng/ml)	CV (%)
1	2.03 ± 0.16	8	1	0.6 ± 0.1	12
2	6.74 ± 0.37	5.6	2	4.6 ± 0.3	6.3

Linearity		Range (ng/ml)	Range (%)	Mean (%)
	Histamine	0.74 - 8.48	85 - 106	100

Recovery		Serial dilution up to	Range (%)	Mean (%)
	Histamine	1:16	92 - 120	103

## 9. References/Literature

- (1) Yagci et al. TCTP/HRF pathway and angiogenesis: A feasible intercourse in chronic lymphocytic leukaemia. *Leukemia Research*, 37:665-670 (2013)
- (2) Coulson et al. Paracetamol (acetaminophen) attenuates in vitro mast cell and peripheral blood mononucleocyte cell histamine release induced by N-acetylcysteine. *Clinical Toxicology*, 48(2):111-114 (2010)
- (3) Rovere et al. Histamine and Selenium in Lung Cancer. *Anticancer Research*, 26: 2937-2942 (2006)

***For updated literature or any other information please contact your local supplier.***

### **Symbols:**

	Storage temperature		Manufacturer		Contains sufficient for <n> tests
	Expiry date		Batch code		For in-vitro diagnostic use only!
	Consult instructions for use		Content		CE labelled
	Caution		Catalogue number		For research use only!

## APPENDIX 1: – For research use only (RUO) –

Protocol for the quantitative determination of Histamine in cell culture samples.

### 1. **Sample collection and storage**

#### **Cell culture**

Samples can be stored at 2 – 8 °C for up to 6 hours. For longer periods ( $\leq 6$  months) the samples should be stored between -20 °C and -80 °C.

### 2. **Test procedure**

Allow all reagents and samples to reach room temperature and mix thoroughly by gentle inversion before use. Duplicate determinations are recommended.

The binding of the antisera and of the enzyme conjugate and the activity of the enzyme are temperature dependent, and the absorption values may vary if a thermostat is not used. The higher the temperature, the higher the absorption values will be. Varying incubation times will have similar influences on the absorbance. The optimal temperature during the Enzyme Immunoassay is between 20 – 25 °C.

*In case of overflow, read the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to 405 nm*

#### 3.1 **Preparation of reagents**

##### **Wash Buffer**

Dilute the 20 ml Wash Buffer Concentrate with water (deionized, distilled, or ultra-pure) to a final volume of 1000 ml.

Storage: 1 month at 2 – 8 °C

##### **Acylation Reagent**

The Acylation Reagent has a freezing point of 18.5 °C. To ensure that the Acylation Reagent is liquid when being used, it must be ensured that the Acylation Reagent has reached room temperature and forms a homogeneous, crystal-free solution.

#### 3.2 **Sample preparation and acylation**

1.	Pipette <b>100 µl</b> of <b>standards, controls</b> and <b>20 µl</b> of <b>cell culture samples</b> into the respective wells of the <b>Reaction Plate</b> .
2.	Add <b>80 µl</b> of <b>Diluent</b> to the wells with <b>cell culture samples</b> .
3.	Add <b>25 µl</b> of <b>Acylation Reagent</b> ( <i>refer to 6.1</i> ) to all wells.
4.	Pipette <b>200 µl</b> of <b>Acylation Buffer</b> into all wells.
5.	Incubate <b>15 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm). <b>Alternatively without shaker:</b> <i>shake the Reaction Plate shortly by hand and incubate <b>15 min at RT</b> (20 – 25 °C).</i>
	<b>Take 25 µl for the subsequent ELISA</b>

#### 3.3 **Histamine ELISA**

1.	Pipette <b>25 µl</b> of the <b>acylated standards, controls and samples</b> into the wells of the <b>Histamine Microtiter Strips</b> .
2.	Pipette <b>100 µl</b> of the <b>Histamine Antiserum</b> into all wells.
3.	Incubate <b>30 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm). <b>Alternatively without shaker:</b> <i>shake the Histamine Microtiter Strips shortly by hand and incubate for <b>40 min at RT</b> (20 – 25 °C).</i>
4.	Discard or aspirate the contents of the wells and wash the plate <b>3 x</b> by adding <b>300 µl</b> of <b>Wash Buffer</b> , <b>discarding</b> the content and <b>blotting dry each time</b> by tapping the inverted plate on absorbent material.
5.	Pipette <b>100 µl</b> of the <b>Enzyme Conjugate</b> into all wells.
6.	Incubate for <b>10 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm). <b>Alternatively without shaker:</b> <i>incubate for <b>20 min at RT</b> (20 – 25 °C).</i>

7. Discard or aspirate the contents of the wells and wash the plate **3 x** by adding **300 µl** of **Wash Buffer**, **discarding** the content and **blotting dry each time** by tapping the inverted plate on absorbent material.
8. Pipette **100 µl** of the **Substrate** into all wells.
9. Incubate for **15 ± 2 min** at **RT** (20 – 25 °C) on a **shaker** (approx. 600 rpm).  
**Alternatively without shaker: incubate for 15 ± 2 min at RT (20 – 25 °C).**
- ⚠ **Avoid exposure to direct sunlight!**
10. Add **100 µl** of the **Stop Solution** to each well and shake the microtiter plate shortly by hand to ensure a homogeneous distribution of the solution.
11. **Read** the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to **450 nm** (if available a reference wavelength between 620 nm and 650 nm is recommended).

#### 4. Calculation of results

Measuring range	Histamine	
	Cell culture	0.6 – 250 ng/ml

The standard curve is obtained by plotting the absorbance readings (calculate the mean absorbance) of the standards (linear, y-axis) against the corresponding standard concentrations (logarithmic, x-axis). Use a non-linear regression for curve fitting (e.g. spline, 4- parameter, akima).

⚠ *This assay is a competitive assay. This means: the OD-values are decreasing with increasing concentrations of the analyte. OD-values found below the standard curve correspond to high concentrations of the analyte in the sample and have to be reported as being positive.*

#### Controls

The concentrations of the **controls** can be read directly from the standard curve.

#### Cell culture samples

The read concentrations of **histamine in cell culture samples** have to be **multiplied by 5**.

The total amount of Histamine excreted in urine during 24 h is calculated as following:

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

## 9. Referenzen/Literatur

- (1) Yagci et al. TCTP/HRF pathway and angiogenesis: A feasible intercourse in chronic lymphocytic leukaemia. Leukemia Research, 37:665-670 (2013)
- (2) Coulson et al. Paracetamol (acetaminophen) attenuates in vitro mast cell and peripheral blood mononucleocyte cell histamine release induced by N-acetylcysteine. Clinical Toxicology, 48(2):111-114 (2010)
- (3) Rovere et al. Histamine and Selenium in Lung Cancer. Anticancer Research, 26: 2937-2942 (2006)

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