IMMUNOASSAYS AND SERVICES

BIOGENIC AMINES & NEUROSCIENCE | ENDOCRINOLOGY | FOOD SAFETY

LABOR DIAGNOSTIKA NORD GmbH & Co.KG | Am Eichenhain 1 | 48531 Nordhorn | Germany | Tel. +49 5921 8197-0 | Fax +49 5921 8197-222 | info@ldn.de | www.ldn.de

Instructions for use
Metanephrine Plasma ELISA **Fast Track**

REF

BA E-8100

IVD

1. Introduction

1.1 Intended use and principle of the test

Enzyme Immunoassay for the quantitative determination of free Metanephrine in plasma.

First, the plasma proteins are removed by precipitation. After this Metanephrine (Metadrenaline) is quantitatively acylated.

The subsequent competitive ELISA uses the microtiter plate format. The antigen is bound to the solid phase of the microtiter plate. The acylated standards, controls and samples and the solid phase bound analytes compete for a fixed number of antibody binding sites. After the system is in equilibrium, free antigen and free antigen-antibody complexes are removed by washing. The antibody bound to the solid phase is detected by an anti-rabbit IgG-peroxidase conjugate using TMB as a substrate. The reaction is monitored at 450 nm.

Quantification of unknown samples is achieved by comparing their absorbance with a reference curve prepared with known standards.

⚠ *The antibodies used in this test kit only recognise the biologically relevant L-forms of Metanephrines. Commercially available synthetic Normetanephrine or Metanephrine is always a mixture of the D- and L-forms. The ratio between both forms differs widely from lot to lot. This has important implications if synthetic Metanephrines are used to enrich native samples. As only about 50% of the synthetic Metanephrines – the L-portion – will be detected by use of this kit, spiked samples will be underestimated. Therefore native samples containing solely the L-form should be used.*

1.2 Clinical application

Metanephrine and Normetanephrine are the metabolites of the catecholamines Epinephrine and Norepinephrine, respectively. Cells derived from neuroendocrine tumors (e.g. pheochromocytoma) are known to produce catecholamines which are secreted episodically via vesicles into the blood stream. But beside this a small portion of the catecholamines is metabolized inside the cells to the corresponding catecholamines metabolites – namely Metanephrine, Normetanephrine and 3-Methoxytyramine – which are secreted at low levels continuously into the blood stream.

Recent studies and publications have shown that the quantification of these plasma free Metanephrine and plasma free Normetanephrine is the most accurate biochemical marker for the clinical diagnosis of pheochromocytoma and follow-up of pheochromocytoma patients.

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as under point "Procedural cautions, guidelines and warnings". Any laboratory result is only a part of the total clinical picture of the patient.

Only in cases where the laboratory results are in an acceptable agreement with the overall clinical picture of the patient it can be used for therapeutic consequences.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

2. Procedural cautions, guidelines, warnings and limitations

2.1 Procedural cautions, guidelines and warnings

- (1) This kit is intended for professional use only. Users should have a thorough understanding of this protocol for the successful use of this kit. Only the test instruction provided with the kit is valid and has to be used to run the assay. Reliable performance will only be attained by strict and careful adherence to the instructions provided.
- (2) This assay was validated for a certain type of sample as indicated in *Intended Use* (please refer to Chapter 1). Any off-label use of this kit is in the responsibility of the user and the manufacturer cannot be held liable.
- (3) Reagents of this kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
- (4) The principles of Good Laboratory Practice (GLP) have to be followed.
- (5) In order to reduce exposure to potentially harmful substances, wear lab coats, disposable protective gloves and protective glasses where necessary.
- (6) All kit reagents and specimens should be brought to room temperature and mixed gently but thoroughly before use. Avoid repeated freezing and thawing of reagents and specimens.
- (7) For dilution or reconstitution purposes, use deionized, distilled, or ultra-pure water.
- (8) The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch with desiccant and used in the frame provided.
- (9) Duplicate determination of sample is highly recommended to be able to identify potential pipetting errors.
- (10) Once the test has been started, all steps should be completed without interruption. Make sure that the required reagents, materials and devices are prepared ready at the appropriate time.
- (11) Incubation times do influence the results. All wells should be handled in the same order and time intervals.

- (12) To avoid cross-contamination of reagents, use new disposable pipette tips for dispensing each reagent, sample, standard and control.
- (13) A standard curve must be established for each run.
- (14) The controls should be included in each run and fall within established confidence limits. The confidence limits are listed in the QC-Report.
- (15) Do not mix kit components with different lot numbers within a test and do not use reagents beyond expiry date as shown on the kit labels.
- (16) Avoid contact with Stop Solution containing 0.25 M H₂SO₄. It may cause skin irritation and burns. In case of contact with eyes or skin, rinse off immediately with water.
- (17) TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them.
- (18) For information on hazardous substances included in the kit please refer to Material Safety Data Sheet (MSDS). The Material Safety Data Sheet for this product is made available directly on the website of the manufacturer or upon request.
- (19) The expected reference values reported in this test instruction are only indicative. It is recommended that each laboratory establishes its own reference intervals.
- (20) The results obtained with this test kit should not be taken as the sole reason for any therapeutic consequence (e.g. medication before a scheduled surgery) but have to be correlated to other diagnostic tests and clinical observations.
- (21) Kit reagents must be regarded as hazardous waste and disposed according to national regulations.

2.2 Limitations

Any inappropriate handling of samples or modification of this test might influence the results.

2.2.1 Interfering substances

Serum/Plasma

Samples containing precipitates or fibrin strands or which are haemolytic or lipemic might cause inaccurate results.

2.2.2 Drug interferences

There are no known substances (drugs) which ingestion interferes with the measurement of Metanephrine level in the sample.

2.2.3 High-Dose-Hook effect

No hook effect was observed in this test.

3. Storage and stability

Store the unopened reagents at 2 - 8 °C until expiration date. Do not use components beyond the expiry date indicated on the kit labels. Once opened the reagents are stable for 1 month when stored at 2 - 8 °C. Once the resealable pouch has been opened, care should be taken to close it tightly with desiccant again.

4. Materials

4.1 Content of the kit

BA D-0023	REAC-TUBES	Reaction Tubes - Ready to use
Content:	Reaction Tubes in a resealable pouch.	
Volume:	2 x 50 tubes	
BA D-0090	FOILS	Adhesive Foil - Ready to use
Content:	Adhesive Foils in a resealable pouch.	
Volume:	1 x 4 foils	
BA E-0030	WASH-CONC 50x	Wash Buffer Concentrate - Concentrated 50x
Content:	Buffer with a non-ionic detergent and physiological pH.	
Volume:	1 x 20 ml/vial, light purple cap	
BA E-0040	CONJUGATE	Enzyme Conjugate - Ready to use
Content:	Goat anti-rabbit immunoglobulins conjugated with peroxidase.	
Volume:	1 x 12 ml/vial, red cap	

Standards and Controls - Ready to use

Cat. no.	Component	Colour/Cap	Concentration	Concentration	Volume/ Vial
			pg/ml	pmol/l	
			MN	MN	
BA R-8301	<u>STANDARD</u> A	white	0	0	12 ml
BA R-8302	<u>STANDARD</u> B	light yellow	36	183	4 ml
BA R-8303	<u>STANDARD</u> C	orange	120	608	4 ml
BA R-8304	<u>STANDARD</u> D	dark blue	360	1 825	4 ml
BA R-8305	<u>STANDARD</u> E	light grey	1 200	6 084	4 ml
BA R-8306	<u>STANDARD</u> F	black	3 600	18 252	4 ml
BA R-8351	<u>CONTROL</u> 1	light green	Refer to QC-Report for expected value and acceptable range!		4 ml
BA R-8352	<u>CONTROL</u> 2	dark red			4 ml

Conversion: Metanephrine (pg/ml) x 5.07 = Metanephrine (pmol/l)

Content: Buffer with stabilizer and a precipitating reagent spiked with a defined quantity of Metanephrine

Hazards identification: 

H302 Harmful if swallowed.

BA E-0055 SUBSTRATE **Substrate** - Ready to use

Content: Chromogenic substrate containing tetramethylbenzidine, substrate buffer and hydrogen peroxide

Volume: 1 x 12 ml/vial, black cap

BA E-0080 STOP-SOLN **Stop Solution** - Ready to use

Content: 0.25 M sulfuric acid.

Volume: 1 x 12 ml/vial, light grey cap

BA E-0131 ADR MN **Metanephrine Microtiter Strips** - Ready to use

Content: 1 x 96 well (12x8) antigen precoated microwell plate in a resealable blue pouch with desiccant

BA E-8110 MN-AS **Metanephrine Antiserum** - Ready to use

Content: Rabbit anti- Metanephrine antibody, blue coloured.

Volume: 1 x 6 ml/vial, blue cap

BA E-8313 ASSAY-BUFF **Assay Buffer** - Ready to use

Content: TRIS-buffer containing proteins and a non-mercury preservative.

Volume: 1 x 12 ml/vial, orange cap

Hazards identification: H317 May cause an allergic skin reaction.

BA R-0028 EQUA-REAG **Equalizing Reagent** - Lyophilized

Content: Human serum, negative for HIV I/II, HBsAg and HCV

Volume: 2 vials, dark green cap

BA R-8312 ACYL-CONC **Acylation Concentrate** - Concentrated

Content: Acylation reagent in DMSO

Volume: 1 x 1.5 ml/vial, dark grey cap

Hazards identification: 

H302 Harmful if swallowed.

4.2 Additional materials and equipment required but not provided in the kit

- Calibrated precision pipettes to dispense volumes between 25 - 500 µl; 3 ml; 10 ml
- Microtiter plate washing device (manual, semi-automated or automated)
- ELISA reader capable of reading absorbance at 450 nm and if possible 620 - 650 nm
- Centrifuge capable of at least 3.000 x g
- microtiter plate shaker (shaking amplitude 3 mm; approx. 600 rpm)
- Absorbent material (paper towel)
- Water (deionized, distilled, or ultra-pure)
- Vortex mixer

5. Sample collection and storage

EDTA-, Citrate- or Heparin-Plasma

Whole blood should be collected into centrifuge tubes (Monovette™ or Vacuette™) containing EDTA, heparin or citrate as anti-coagulant and centrifuged (according to manufacturer's instructions) immediately after collection.

Haemolytic and lipemic samples should not be used for the assay.

Storage: up to 6 hours at 2 - 8 °C, for longer period (up to 6 month) at -20 °C.

Repeated freezing and thawing should be avoided.

6. Test procedure

Allow all reagents to reach room temperature and mix thoroughly by gentle inversion before use. Number the Reaction Tubes accordingly. Duplicate determinations are recommended.

The binding of the antibodies and the enzyme conjugates and the activity of the enzyme used are temperature dependent, and the absorption values may vary if a thermostat is not used. The higher the temperature, the higher the absorption values will be. The absorption values also depend on the incubation times. The optimal temperature for the Enzyme Immunoassay is between 20 - 25 °C.

⚠ *In case of overflow, read the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to 405 nm.*

6.1 Preparation of reagents

Wash Buffer

Dilute the 20 ml of Wash Buffer Concentrate with water (deionized, distilled, or ultra-pure) to a final volume of 1000 ml.

Storage: 1 month at 2 - 8 °C

Equalizing Reagent

Reconstitute the Equalizing Reagent with 10 ml water (deionized, distilled, or ultra-pure).

Reconstituted Equalizing Reagent which is not used immediately has to be stored in aliquots for max. 1 month at -20 °C and may be thawed only once.

Acylation Solution

As the Acylation Solution is only stable for a maximum of 3 minutes it should not be prepared before starting the assay. Therefore its preparation is described in the protocol in chapter 6.3, step 3.

Discard after use!

6.2 Precipitation

1. Pipette **100 µl** of **standards**, **100 µl** of **controls**, and **500 µl** of **plasma samples** into the respective **Reaction Tubes**.

2. Add **500 µl Equalizing Reagent** (*refer to 6.1*) to all tubes containing standards and controls.

3. Add **100 µl Standard A** to all tubes containing plasma samples.

4. Mix **Reaction Tubes** thoroughly (vortex) and centrifuge for **15 minutes** at **3,000 x g**.

⚠ Take **75 µl** of the clear supernatant for the Metanephrine ELISA.

6.3 Metanephrine ELISA

1.	Pipette 50 µl of Assay Buffer into the appropriate wells of the Metanephrine Microtiter Strips .
2.	Pipette 75 µl of the clear supernatant from the standards, controls and samples into the wells.
3.	Preparation of Acylation Solution : Pipette 80 µl Acylation Reagent Concentrate (BA R-8312) to 3 ml water (deionized, distilled, or ultra-pure) and mix thoroughly.
4.	Pipette 25 µl of the freshly prepared Acylation Solution into all wells.
5.	Incubate for 15 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm).
6.	Pipette 50 µl of the Metanephrine Antiserum into all wells.
7.	Cover the plate with Adhesive Foil . Shake for 1 min at RT (20 – 25 °C) on a shaker and incubate for 15 - 20 h (overnight) at 2 – 8 °C without shaking.
8.	Remove the foil. Discard or aspirate the content of the wells. Wash the plate 4 x by adding 300 µl of Wash Buffer , discarding the content and blotting dry each time by tapping the inverted plate on absorbent material.
9.	Pipette 100 µl of the Enzyme Conjugate into all wells.
10.	Incubate for 30 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm).
11.	Discard or aspirate the content of the wells and wash the plate 4 x by adding 300 µl of Wash Buffer , discarding the content and blotting dry each time by tapping the inverted plate on absorbent material.
12.	Pipette 100 µl of the Substrate into all wells and incubate for 20 - 30 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). Avoid exposure to direct sunlight!
13.	Add 100 µl of the Stop Solution to each well and shake the microtiter plate to ensure a homogeneous distribution of the solution.
14.	Read the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to 450 nm (if available a reference wavelength between 620 nm and 650 nm is recommended).

7. Calculation of results

Measuring range	Metanephrine
	17 – 3 600 pg/ml

The standard curve is obtained by plotting the absorbance readings (calculate the mean absorbance) of the standards (linear, y-axis) against the corresponding standard concentrations (logarithmic, x-axis). Use a non-linear regression for curve fitting (e.g. spline, 4- parameter, akima).

 *This assay is a competitive assay. This means: the OD-values are decreasing with increasing concentrations of the analyte. OD-values found below the standard curve correspond to high concentrations of the analyte in the sample and have to be reported as being positive.*

The concentrations of the **samples** and **controls** can be read directly from the standard curve.

Samples found with concentrations higher than the highest standard (Standard F) should be diluted accordingly with Equalizing Reagent (BA R-0028) and have to be re-assayed.

Conversion

Metanephrine (pg/ml) x 5.07 = Metanephrine (pmol/l)

Expected reference value

It is strongly recommended that each laboratory should determine its own reference value.

Metanephrine
< 90 pg/ml


7.1 Quality control

It is recommended to use control samples according to national regulations. Use controls at both normal and pathological levels. The kit or other commercial controls should fall within established confidence limits. The confidence limits of the kit controls are indicated in the QC-Report.

7.2 Typical standard curve

Example, do not use for calculation!

8. Assay characteristics

Analytical Sensitivity (Limit of Detection)		Metanephrine
	Plasma	17 pg/ml

Analytical Specificity (Cross Reactivity)	Substance	Cross Reactivity (%)
		Metanephrine
	Derivatized Metanephrine	100
	Derivatized Normetanephrine	0.04
	3-Methoxytyramin.HCl	< 0.001
	Adrenaline	< 0.001
	Noradrenaline	< 0.001
	Dopamin.HCl	< 0.001
	VMS	< 0.001
	HMVS	< 0.001
	L-DOPA	< 0.001
	L-Tyrosin	< 0.001
	Tyramine.HCl	< 0.001
	Normetanephrine	< 0.001
	Acetaminophen	< 0.001

Precision							
Intra-Assay				Inter-Assay			
	Sample	Range (pg/ml)	CV (%)		Sample	Range (pg/ml)	CV (%)
Metanephrine	1	155 ± 17	11	Metanephrine	1	133 ± 13	10
	2	245 ± 28	11		2	257 ± 23	8.9
	3	523 ± 38	7.3		3	528 ± 50	9.5

Linearity			Range	Serial dilution up to	Mean (%)
	Metanephrine	Plasma	40 – 2100	1: 65	99

Recovery			Mean (%)	Range (%)	% Recovery after spiking
	Metanephrine	Plasma	97	85 – 111	


Method Comparison ELISA vs. LC-MS/MS	Metanephrine	Plasma	LC-MS/MS = 1.1x - 19.4	r = 0.98; n = 59
---	--------------	--------	------------------------	------------------

9. References/Literature

- (1) Jeyaraman et al. The role of urinary fractionated metanephrines in the diagnosis of pheochromocytoma. The Indian Journal of Medical Research, 137(2):316-323 (2013)
- (2) Unger et al. Plasma and Urinary Metanephrines Determined by an Enzyme Immunoassay, but not Serum Chromogranin A for the Diagnosis of Pheochromocytoma in Patients with Adrenal Mass. Exp Clin Endocrinol Diabetes, 120:494-500 (2012)
- (3) Stefanescu et al. Salivary free catecholamines metabolites as possible biochemical markers in pheochromocytoma diagnosis. Acta Endocrinologica (Buc), VII (4): 431-442 (2011)

 **For updated literature or any other information please contact your local supplier.**

Symbols:

	Storage temperature		Manufacturer		Contains sufficient for <n> tests
	Expiry date	LOT	Batch code	IVD	For in-vitro diagnostic use only!
	Consult instructions for use	CONT	Content	CE	CE labelled
	Caution	REF	Catalogue number	RUO	For research use only!

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