

Instructions for use
DHEA Saliva ELISA

REF

SA E-7000R



96

RUO

For research
use only –
Not for use
in diagnostic
procedures

DHEA Saliva ELISA

1. INTRODUCTION

1.1 Intended use and principle of the test

An enzyme immunoassay for the measurement of active free DHEA in saliva.

2. PRINCIPLE

The *DHEA Saliva ELISA* Kit is a solid phase enzyme-linked immunosorbent assay (ELISA), based on the principle of competitive binding. The microtiter wells are coated with an antibody directed against the DHEA molecule. Endogenous DHEA of a patient sample competes with a DHEA-horseradish peroxidase conjugate for binding to the coated antibody. After incubation the unbound conjugate is washed off.

The amount of bound peroxidase conjugate is inversely proportional to the concentration of DHEA in the sample. After addition of the substrate solution, the intensity of color developed is inversely proportional to the concentration of DHEA in the patient sample.

3. WARNINGS and PRECAUTIONS

1. This kit is for research use only.
2. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.
3. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
4. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
5. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
6. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
7. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
8. Allow the reagents to reach room temperature (21 - 26 °C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the patient samples will not be affected.
9. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
10. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
11. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
12. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
13. Do not use reagents beyond expiry date as shown on the kit labels.
14. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiterplate readers.
15. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
16. Avoid contact with Stop Solution. It may cause skin irritation and burns.
17. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
18. For information please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from the manufacturer.

4. REAGENTS

4.1 Reagents provided

SA E-7031	 96	Microtiterwells
Content:	Wells coated with an anti-DHEA antiserum (polyclonal).	
Volume:	12 x 8 (break apart) strips, 96 wells.	

Standards and Controls -Ready to use

Cat.-no.	Component	Concentration	Volume / Vial
SA E-7001	STANDARD A	0 pg/ml	3.0 ml
SA E-7002	STANDARD B	10 pg/ml	1.0 ml
SA E-7003	STANDARD C	40 pg/ml	1.0 ml
SA E-7004	STANDARD D	160 pg/ml	1.0 ml
SA E-7005	STANDARD E	640 pg/ml	1.0 ml
SA E-7006	STANDARD F	2560 pg/ml	1.0 ml
SA E-7051	CONTROL 1	For control values and ranges please refer to QC-Datasheet.	1.0 ml
SA E-7052	CONTROL 2		1.0 ml

SA E-7040 CONJUGATE Enzyme Conjugate - Ready to use

Content: DHEA conjugated to horseradish peroxidase.

Volume: 1 x 11 ml/vial

AR E-0055 SUBSTRATE Substrate Solution - Ready to use

Content: Tetramethylbenzidine (TMB).

Volume: 1 x 22 ml/vial

Hazards identification:



H360D May damage the unborn child.

AR E-0080 STOP-SOLN Stop Solution - Ready to use

Content: contains 2 N hydrochlorid acid solution.

Volume: 1 x 7 ml/vial

avoid contact with the stop solution. It may cause skin irritations and burns.

Hazards identification:



H290 May be corrosive to metals.

H314 Causes severe skin burns and eye damage.

H335 May cause respiratory irritation.

AR E-0030 WASH-CONC 10x Wash Solution - Concentrated 10x

Volume: 1 x 50 ml/vial

see "Preparation of Reagents"

All reagents contain azide-free and mercury-free preservatives.

Note: Additional Standard A for sample dilution is available upon request.

4.2 Material required but not provided

- microcentrifuge
- A microtiterplate reader capable for endpoint measurement at 450 ± 10 nm
- Microplate mixer operating at about 600 - 900 rpm
- Vortex mixer
- Calibrated variable precision micropipettes (50 μ l, 100 μ l, 200 μ l)
- Absorbent paper
- Distilled or deionized water
- Timer
- Semi logarithmic graph paper or software for data calculation

4.3 Storage conditions

When stored at 2°C to 8°C unopened reagents will be stable until expiration date. Do not use reagents beyond this date. Opened reagents must be stored at 2° - 8 °C. After first opening the reagents are stable for 30 days if used and stored properly.

Microtiter wells must be stored at 2°C to 8°C. Take care that the foil bag is sealed tightly.

4.4 Reagents preparation

Allow the reagents and the required number of wells to reach room temperature (21 - 26 °C) before starting the test.

Wash Solution

Add deionized water to the 10X concentrated Wash Solution.

Dilute 50 ml of concentrated Wash Solution with 450 ml deionized water to a final volume of 500 ml.

The diluted Wash Solution is stable for 3 months at room temperature.

4.5 Disposal of the kits

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Safety Data Sheet.

4.6 Damaged test kits

In case of any severe damage of the test kit or components, the manufacturer has to be informed written, latest one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

5. SPECIMEN COLLECTION AND PREPARATION

Samples containing sodium azide should not be used in the assay. The saliva samples should be completely colorless. Even the slightest red color shows blood contamination. Such blood contamination will give falsely elevated concentration values. In case of visible blood contamination the patient should discard the sample, rinse the collection device with water, also rinse the mouth with (preferably) cold water, wait for 10 minutes and take a new sample. Do not chew anything during the sampling period. Any pressure on the teeth may result in falsely elevated measurements due to an elevated content of gingival liquid in the saliva sample.

5.1 Specimen Collection

For the correct collection of saliva we are recommending to use only appropriate devices made from ultra-pure polypropylene. Do not use any PE devices or cotton based Salivettes for sampling. False readings will result. Glass tubes can be used as well, but in this case special attention is necessary for excluding any interference caused by the stopper. Please contact the manufacturer for more details.

As food might contain significant amounts of steroid hormones samples preferably should be taken while fasting. If fasting should be a problem at least any food of animal origin (meat or dairy products) should be avoided prior to finalizing the collection. In the morning breakfast should be done only after finalizing the collection procedure. During the day the collection period should be timed just before an anticipated meal. As the steroid hormone secretion in saliva as well as in serum shows an obvious dynamic secretion pattern throughout the day it is important to always collect 5 samples during a 2 hour period; this means every 30 minutes one sample. If possible the volume of each single sample should be a minimum of 0.5 ml (better 1 ml). Saliva flow may be stimulated by drinking water. This is allowed and even recommended before and during the collection period. Drinking of water is not allowed during the last 5 minutes before taking the single samples.

5.2 Specimen Storage and Preparation

Saliva samples in general are stable at ambient temperature for several days. Therefore mailing of such samples by ordinary mail without cooling will not create a problem. Storage at 4 °C can be done for a period of up to one month. Whenever possible samples preferable should be kept at a temperature of -20 °C. Even repeated thawing and freezing is no problem. Each sample has to be frozen, thawed, and centrifuged at least once anyhow in order to separate the mucins by centrifugation. Upon arrival of the samples in the lab the samples have to be kept frozen at least overnight. Next morning the frozen samples are warmed up to room temperature and mixed carefully. Then the samples have to be centrifuged for 5 to 10 minutes. Now the clear colorless supernatant is easy to pipette. If the sample should show even a slightly red color it should be discarded. Otherwise the concentration value most probably will be falsely elevated. Due to the episodic variations of the steroid secretion we highly recommend the strategy of multiple sampling. If such a set of multiple samples has to be tested the lab (after at least one freezing, thawing, and centrifugation cycle) has to mix the aliquots of the 5 single samples in a separate sampling device and perform the testing from this mixture.

5.3 Specimen Dilution

If in an initial assay, a specimen is found to contain more than the highest standard, the specimens can be diluted with Standard A and re-assayed as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

- a) Dilution 1:10: 10 µl saliva + 90 µl Standard A (mix thoroughly)
- b) Dilution 1:100: 10 µl of dilution a) + 90 µl Standard A (mix thoroughly).

6. ASSAY PROCEDURE

6.1 General remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination.
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.
- Respect the incubation times as stated in this instructions for use

6.2 Assay Procedure

Each run must include a standard curve.

1.	Secure the desired number of coated strips in the frame holder.
2.	Dispense 100 µl of each Standard, Control and sample <u>with new disposable tips</u> into appropriate wells.
3.	Dispense 100 µl of Enzyme Conjugate into each well.
4.	Incubate for 60 minutes at room temperature on a Microplate mixer (≥ 600 rpm). <i>Important note: Optimal reaction in this assay is markedly dependent on shaking of the microplate!</i>
5.	Briskly empty the contents of the wells by aspiration or by decanting. Rinse the wells 4 times with diluted Wash Solution (300 µl per well). Strike the wells sharply on absorbent paper to remove residual droplets. <i>Important note: The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!</i>
6.	Add 200 µl of Substrate Solution to each well.
7.	Incubate for 30 minutes in the dark at room temperature.
8.	Stop the enzymatic reaction by adding 50 µl of Stop Solution to each well.
9.	Determine the absorbance of each well at 450±10 nm . It is recommended that the wells are read <u>within 15 minutes</u> .

6.3 Calculation of results

1. Calculate the average absorbance values for each set of standards, controls and patient samples.
2. Construct a standard curve by plotting the mean absorbance obtained from each standard against its concentration with absorbance value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the package insert have been calculated automatically using a 4 PL (4 Parameter Logistics) curve fit. 4 Parameter Logistics is the preferred calculation method. Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this standard curve. Samples with concentrations higher than that of the highest standard have to be further diluted. For the calculation of the concentrations this dilution factor has to be taken into account.

6.3.1 Example of typical Standard Curve

The following data is for demonstration only and cannot be used in place of data generations at the time of assay.

Standard	Optical Units (450 nm)
Standard A 0 pg/ml	2.940
Standard B 10 pg/ml	2.701
Standard C 40 pg/ml	2.290
Standard D 160 pg/ml	1.657
Standard E 640 pg/ml	0.890
Standard F 2560 pg/ml	0.426

7. EXPECTED NORMAL VALUES

In order to determine the normal range of DHEA free in Saliva samples from adult male and female apparently healthy subjects were collected and analyzed using the *DHEA Saliva ELISA* kit. The following range was calculated from this study.

Age Group (years)	Men			Women		
	5% - 95% Percentile (pg/ml)	Median (pg/ml)	n	5% - 95% Percentile (pg/ml)	Median (pg/ml)	n
< 21	30.4 - 537.7	200.7	7	27.2 - 565.4	215.7	24
21 - 30	291.4 - 826.7	464.4	10	73.5 - 780.7	605.2	50
31 - 40	306.7 - 892.3	514.2	10	124.5 - 745.1	335.0	50
41 - 50	86.8 - 713.7	285.2	25	85.7 - 480.8	222.3	50
51 - 60	79.1 - 525.3	228.4	23	76.7 - 620.2	217.7	50
> 60	39.4 - 694.9	171.2	28	34.7 - 467.1	170.8	50

8. QUALITY CONTROL

Good laboratory practice requires that controls need to be run with each Standard curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or the manufacturer directly.

9. PERFORMANCE CHARACTERISTICS

9.1 Analytical Sensitivity

The analytical sensitivity of the *DHEA Saliva ELISA* was calculated by subtracting 2 standard deviations from the mean of twenty-two (22) replicate analyses of Standard A.

The analytical sensitivity of the assay is 3.7 pg/ml.

9.2 Specificity (Cross Reactivity)

The following materials have been evaluated for cross reactivity.

Steroids	% Crossreactivity at 50% Binding
DHEA-S	< 0.01
Testosterone	< 0.01
5 α -Dihydrotestosterone	< 0.01
Androstendione	0.06
Progesterone	0.23
17 α -Hydroxyprogesterone	< 0.01
Pregnenolone	0.01
17-Hydroxy-Pregnenolone	0.07
Deoxycorticosterone	0.05
Corticosterone	< 0.01
Cortisol	< 0.01
11-Desoxycortisol	0.01
Estradiol-17 β	< 0.01
Estradiol-17 α	< 0.01
Estrone	< 0.01
Estriol	< 0.01

9.3 Assay Dynamic Range

The range of the assay is between 10 - 2560 pg/ml.

9.4 Reproducibility

9.4.1 Intra-Assay

The intra-assay variation was determined by replicate measurements of 3 saliva samples within one run using the *DHEA Saliva ELISA*. The within-assay variation is shown below:

	Sample 1	Sample 2	Sample 3
Mean (pg/ml)	255.1	672.8	822.8
SD (pg/ml)	13.3	60.7	39.9
CV (%)	5.2	9.0	4.9
n =	19	19	19

9.4.2 Inter-Assay

The inter-assay variation was determined by duplicate measurements of 3 saliva samples in 10 different runs using the *DHEA Saliva ELISA*. The inter-assay variation is shown below:

	Sample 1	Sample 2	Sample 3
Mean (pg/ml)	238.6	648.8	797.5
SD (pg/ml)	10.6	44.2	51.2
CV (%)	4.5	6.8	6.4
n =	10	10	10

9.5 Recovery

Using the standard matrix a spiking solution of 20 ng DHEA/ml was prepared. 500 µl of three saliva were spiked with 5, 10 and 15 µl of the spiking solution leaving the saliva matrices relatively intact. All samples were measured by the *DHEA Saliva ELISA* procedure.

Saliva	Spiking (pg/ml)	Observed (pg/ml)	Expected (pg/ml)	Recovery (%)
1	-	291	-	-
	200	517	491	105
	400	668	691	97
	600	859	891	96
2	-	64	-	-
	200	287	264	109
	400	520	464	112
	600	661	664	100
3	-	247	-	-
	200	490	447	110
	400	668	647	103
	600	847	847	100

9.6 Linearity

Four saliva samples containing different amounts of analyte were serially diluted with Standard A and assayed with the *DHEA Saliva ELISA*. The percentage recovery was calculated by comparing the expected and observed values for DHEA.

Saliva	Dilution	Observed (pg/ml)	Expected (pg/ml)	Linearity (%)
1	native	569	-	-
	1 in 2	284	285	100
	1 in 4	127	142	89
	1 in 8	77	71	108
2	native	384	-	-
	1 in 2	167	192	87
	1 in 4	78	96	81
	1 in 8	40	48	83
3	native	238	-	-
	1 in 2	86	119	72
	1 in 4	63	60	105
	1 in 8	30	30	100
4	native	292	-	-
	1 in 2	135	146	92
	1 in 4	73	73	100
	1 in 8	40	37	108

10. LIMITATIONS OF PROCEDURE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice. Any improper handling of samples or modification of this test might influence the results.

10.1 High-Dose-Hook Effect

No hook effect was observed in this test.

10.2 Drug Interference

Any medication (cream, oil, pill etc) containing DHEA of course will significantly influence the measurement of this analyte in saliva.

11. REFERENCES

1. Assessing cortisol and dehydroepiandrosterone (DHEA) in saliva: effects of collection method.
Gallagher P, Leitch MM, Massey AE, McAllister-Williams RH, Young AH
J Psychopharmacol, Sep 2006 (Vol. 20, Issue 5, Pages 643-9)
2. Effects of DHEA administration on episodic memory, cortisol and mood in healthy young men: a double-blind, placebo-controlled study.
Alhaj HA, Massey AE, McAllister-Williams RH Psychopharmacology (Berl), Nov 2006 (Vol. 188, Issue 4, Pages 541-51)
3. Bacteria in the oral mucosa and its effects on the measurement of cortisol, dehydroepiandrosterone, and testosterone in saliva.
Whembolua GL, Granger DA, Singer S, Kivlighan KT, Marguin JA Horm Behav, Apr 2006 (Vol. 49, Issue 4, Pages 478-83)
4. Anthropometry and body composition do not predict bioavailable androgen or progesterone concentration in adolescent girls.
Bond LJ, Vella ET, Kiparissis Y, Wynne-Edwards KE Am J Hum Biol, Sep 2006 (Vol. 18, Issue 5, Pages 639-53)
5. The influence of 10 min of the Johrei healing method on laboratory stress.
Laidlaw TM, Naito A, Dwivedi P, Hansi NK, Henderson DC, Gruzelier JH Complement Ther Med, Jun 2006 (Vol. 14, Issue 2, Pages 127-32)
6. Salivary cortisol, dehydroepiandrosterone-sulphate (DHEA-S) and testosterone in women with chronic migraine.
Patacchioli FR, Monnazzi P, Simeoni S, De Filippis S, Salvatori E, Coloproisco G, Martelletti P J Headache Pain, Apr 2006 (Vol. 7, Issue 2, Pages 90-4)
7. Aggression, dominance, and affiliation: Their relationships with androgen levels and intelligence in 5-year-old children.
Azurmendi A, Braza F, Garcia A, Braza P, Munoz JM, Sanchez-Martin JR Horm Behav, Jun 2006 (Vol. 50, Issue 1, Pages 132-40)
8. Cognitive abilities, androgen levels, and body mass index in 5-year-old children.
Azurmendi A, Braza F, Sorozabal A, Garcia A, Braza P, Carreras MR, Munoz JM, Cardas J, Sanchez-Martin JR Horm Behav, Aug 2005 (Vol. 48, Issue 2, Pages 187-95)
9. Effect of prolonged stress on the adrenal hormones of individuals with irritable bowel syndrome.
Sugaya N., Izawa S., Saito K., Shiotsuki K., Nomura S. and Shimada H.
BioPsychoSocial Medicine 2015, 9:4

⚠ For updated literature or any other information please contact your local supplier.

Symbols:

	Storage temperature		Manufacturer		Contains sufficient for <n> tests
	Expiry date	LOT	Batch code		
	Consult instructions for use	CONT	Content		
	Caution	REF	Catalogue number	RUO	For research use only!